



Electric Multifunctional Operating Table
Fury 10



Fury 10

Features

Lowest Position of
Table Top

710_{mm}



Highest Position of
Table Top

1010_{mm}



Back Plate
Upward

0°-80°



Back Plate Downward

0°-20°



Reversed
Trendelenburg

0°-25°



Trendelenburg

0°-25°



Head Plate Upward

0°-40°



Head Plate Downward

0°-90°



Right Lateral Tilt

0-20°



Left Lateral Tilt

0-20°



Horizontal Sliding

350_{mm}



Leg Plate Rest
-90°~15°




Back Plate Downward

0°-20°



Head Plate Downward

0°-90°



Head Plate Downward

0°-90°



Fury 10



Armrest



Anesthesia Frame



Leg Crutches



Shoulder Support



Body Supporter



Fury 10

Technical Chart

Technical Configuration

Bed Frame	1pc
Mattress	1set
Anesthesia Screen Frame	1pc
Shoulder Holder with Cushion	2pc
Arm Holder with Cushion	2pc
Leg Support with Cushion	2pc
Foot Support with Cushion	1pc
Kidney Bridge Handle	1pc

Electronics

Battery Standby Time(h)	80
Power Supply	220V 50HZ 350VA
Power Control Circuit	12V
Plug	As required

Technical Parameters

External L*W*H	2060*500*(710-1010)mm
Reverse Trendelenburg	0-25°
Trendelenburg	0-25°
Lateral Tilt Left/Right	0-15°
Head Plate Up/Down	0-40° / 0-90°
Back Plate Up/Down	0-80° / 0-20°
Leg Plate Up/Dpwn	0-15°/0-90°
Kidney Bridge Up	110mm
Horizontal Sliding C type Arm	350mm

Quality Control

MATERIAL

The base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned



POWDER SUPPLY

Internal battery maintaining power supply for 50-80 operations and with low battery warning function
50-80




REMOTE PANEL & COLUMN PANEL

Leg Plate With Imported Gas Spring For Easy Control.



BRAKING SYSTEM

Pedal Brake, Easy To Move And Lock.



WHEEL

6" castors import from Germany



HYDRAULIC SYSTEM

Brand hydraulic system (German Hoerbiger) make table with low noise, high safety and reliability



Electro Multifunctional
Operating Table

Fury 10





3403 NW 82nd Avenue, Suite 102 E.
Doral, Florida 33122. USA
www.hosmed-inc.com