



Electric Multifunctional Operating Table
Fury 5



Fury 5

Features

Lowest Position of
Table Top

750mm



Highest Position of
Table Top

1050mm



(Optional) Head Rest



Reversed Trendelenburg

0°-22°



Trendelenburg

0°-25°



(Optional) Lateral Arm
Support



Right Lateral Tilt

0°-20°



Left Lateral Tilt

0°-20°



(Optional) Waste Basin
For Gynecology



Anorectum
Operation Position



E.N.T Operation Position



(Optional) Orthopedics
Traction Frame



Fury 5



Armrest



Anesthesia Frame



I.V. Pole



Leg Crutches



Body Supporter



Fury 5

Technical Chart

Technical Parameters

Reverse Trendelenburg/Trendelenburg	22°/25°
Lateral Tilt Right/Left	20°/20°
Head Rest (Up/Down)	45°/90°
Foot Plate (Up/Down)	15°/90°
Foot Plate Horizontal Positioning	90°
Back Plate (Up/Down)	75°/20°
Kidney Bridge Elevation	110mm
Horizontal Sliding	300mm

Dimension Parameter

External Length	2100mm
External Width	500mm
Table Height Min.	800mm
Table Height Max.	1000mm

Technical Configuration

Bed Frame	1pc
Mattress	1set
Anesthesia Screen Frame	1pc
Shoulder Holder with Cushion	2set
Arm Holder with Cushion	2set
Leg Support with Cushion	2set
Foot Support with Cushion	2set
Kidney Bridge Handle	1pc

Quality Control



Electric Multifunctional
Operating Table

Fury 5





3403 NW 82nd Avenue, Suite 102 E.
Doral, Florida 33122. USA
www.hosmed-inc.com