



# Manual Reclining Dialysis Chair **NDC1**



# Features and Highlights

## Technical parameters

• Seat Size (LxWxH)	820×850×1070 mm
• Bed Size (LxWxH)	820×1600×450 mm

**Design**  
Through the button to achieve a key flat, seconds into a recliner, can lie down to sleep.



**Leather**  
Super strength PU material, Environmental protection.



**Backrest**  
Sectional type backrest gives the patient strong support and effectively protects the patient's spine.



**Sponge**  
Flame retardant, 4 cm thickness, high-density foam,density 30 kg/ m³.



**Armrest**  
Solid wood armrest with rounded edges and a comfortable feel.



**Rotary Table**  
Solid wood rotary table, smooth and comfortable surface.



**Cup Base**  
It is convenient for the patient to put the cup.



**Castors**  
High quality medical silent casters, 3" casters , anti-static electricity, protect the ground, waterproof and anti-bacteria.





# Manual Reclining Dialysis Chair NDC1





3403 NW 82nd Avenue, Suite 102 E.  
Doral, Florida 33122. USA  
[www.hosmed-inc.com](http://www.hosmed-inc.com)