



Electric Dialysis Chair
NDC2



Features

Rich Parts

Leading brand of Medical Furniture



Optional genuine leather or supebric cover.with art1sanal making.comfortable hand feeling



Ergonomics design removable pollow



Environmental plastic can be sterilized 1n l 00°C, tensile upto 30MPa
Adjustable ABS tray table, with ant1sk1d convex edge, tank for cup and cellphone . USB/SV optional.
Very convenient for donor.



Full SS IV pole, bearing 15kg

Technical Chart

QUALITY CONTROL

DESIGN
European design. 4 sections ergonomics design



STEEL
Baosteel from Fortune Global 500



PAINTING
11 Process epoxy painting. ASTM anti-bacterial. paint thickness 0.12mm. brightness 60°. paint can resist 50kg impact
0.12_{mm} | 60° | 50_{kg}





WELDING
Panasonic Robotic ensure 360° full smooth wedding
360°



ELECTRIC SYSTEM
American motor, with UL/ROHS/EN standard. Noise under 50dB, with lifetime 20,000times movement.
50_{dB} | 20,000_{times}



ARMREST
Armrest updated. can be universal adjust. Donor can adjust any position he/she wants
Armrest is R 120° cone avoid donor's arm sliding the donor be safety and comfortable R120°

Back Section Leg Section 560_{mm}
the lowest height is 560mm, easy to sit CPR

MAIN TECHNICAL PARAMETERS

| | |
|--------------------|-----------------|
| Size | 1830mm x 580mm |
| Back | 800*580mm |
| Seat | 530*580mm |
| Leg | 490*580mm |
| Foot board | (230-360)*220mm |
| Seat height | 560mm |
| Pillow | 370*290mm |
| Armrest | S00x 160mm |
| Mattress thickness | 60mm |

FUNCTION

| | |
|-------------------|-----------|
| Back rest lifting | 0° -65° |
| Leg rest lifting | 10° -90° |
| Armrest swing out | 530*580mm |
| With CPR | |

TECHNICAL CONFIGURATION

| | |
|--------------------|-------|
| Hand Controller | 1pc |
| Power Line | 1pc |
| 3" Silent Caster | 1 set |
| Small dining table | 1pc |

ELECTRONICS

| | |
|-----------|------------------------------|
| Voltage | 100-240V |
| Frequency | S0/60Hz |
| Plug | as required (hospital grade) |



3403 NW 82nd Avenue, Suite 102 E.
Doral, Florida 33122. USA
www.hosmed-inc.com