



Electric Dialysis Chair  
**NDC3**



# Features

## Rich Parts

Leading brand of Medical Furniture



Foot pedal connects with leg rest, adjustable, can automatic extension the chair length 250mm, realize both foot step and leg support functions.



Can choose genuine leather or supebnc cover, with art1sanal making, comfortable hand feeling.



Environmental plastic can be sterilized 1n l 00° (, tensile upto 30MPa Adjustable ABS tray table, with ant1sk1d convex edge, tank for cup and cellphone. USB/SV optional. Very convenient for donor.



Environmental plastic can be sterilized 1n l 00° (, tensile upto 30MPa Adjustable ABS tray table, with ant1sk1d convex edge, tank for cup and cellphone. USB/SV optional. Very convenient for donor.



# Technical Chart



Back Section

Leg Section

Height Adjustment



The lowest height is 550mm. easy to sit



CPR

## MAIN TECHNICAL PARAMETERS

Seat Size	1800mm x 580mm (L*W)
Seat Height Adjustable	550mm-750mm
Back Size	770*580mm(L *W)
Seat Size	470*580mm(L*W)
Leg Size	380*540mm(L *W)
Foot Board Size	410*200mm(L *W)
Seat Height	550-750mm
Mattress Thickness	6cm
Back Rest Lifting	0° -70°
Leg Rest Lifting	20°-75°
Armrest Swing Out	0° -270°
Armrest Size	500 x 160mm(L*W)
Armrest Height Adjustable	100mm
Four 3" Castors With Break	

## MAIN TECHNICAL PARAMETERS

Hand Controller	1pc
Power Line	1pc
3" Silent Caster	1pc
Small dining table	1pc
With CPR	

## ELECTRONICS

Voltage	100-240V
Frequency	50/60Hz
3" Silent Caster	as required (hospital grade)

# QUALITY CONTROL

## DESIGN

European design.  
4 sections ergonomics design



## STEEL

Baosteel from Fortune Global 500

## PAINTING

11 Process epoxy painting. ASTM testing anti-baternal. paint thickness 0.12mm. brightness 60°. paint can resist 50kg impact

0.12<sub>mm</sub> | 60° | 50<sub>kg</sub>

WELDING  
PANASONIC  
REBOTIC ENSURE  
360° full smooth  
wedding

360°


50<sub>dB</sub> | 20,000<sub>times</sub>

## ELECTRIC SYSTEM

American motor, with UL/ROHS/EN standard. Noise under 50dB, with lifetime 20,000 times movement.

Armrest updated. can be universal adjust. Donor can adjust any position he/she wants.

Armrest is RI 20° concave curve. avoid donor's arm sliding. ensure the donor be safety and comfortable





3403 NW 82nd Avenue, Suite 102 E.  
Doral, Florida 33122. USA  
[www.hosmed-inc.com](http://www.hosmed-inc.com)