



# Hydraulic Stretcher **NS5H**



# Features



Hydraulic treadle for rise



(1 Leg section lowers down to make head higher (Than Leg Section)  
 2. Back section lowers down to make leg section higher (Than head Section)  
 3. Step on the middle of the pedal, make the bed base fall down steadily



Four double side casters of 6 inch diameter, simultaneously lockable by the step pedal at both sides.

TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts



Waterproof PU mattress, convenient for transfer patient



Foldable aluminium alloy guardrail



Bumper wheel detachable, max. protect the bed and patient during transportation



Full SS IV pole, bearing 15kg

# Technical Chart

## QUALITY CONTROL





## Technical parameters

Size	2100*850*(630-930)mm
Back-rest	0-75° (±10°)
Knee-rest	0-45°
Reverse Trendelenburg	0-15°
Trendelenburg	0-15°

## Technical configuration

USA Hydraulic Pump	2 pcs.
Crank System	1 set
Back Rest Air Spring	2 pcs.
8" German Castors, Central Locking	4 pcs.
Emergency Push Handle	1 set
Folding Side Rail	1 set
74mm Wall Bumper Wheel	1 set
IV pole prevision	2 pcs.
Mattress	1 pcs.
IV Pole	1 set
Central Locking Break Pedal	1 pcs.
X-Ray Translucent Bed Platform	1 set







3403 NW 82nd Avenue, Suite 102 E.  
Doral, Florida 33122. USA  
[www.hosmed-inc.com](http://www.hosmed-inc.com)