




Electric Hydraulic Operating Table  
**Pisa 3**



# Pisa 3

## Features



### Material

The base, column and side hoder are all made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.




### Powder Supply

Internal battery maintaining power supply for 50-80 operations and with low battery warning function.




### Remote Panel

All function controlled by hand controller in the column, and can be locked.




### Hydraulic System

Rand hydraulic system (German Hoerbiger) make table with low noise, high safety and reliability.



### Control System

Foot treadle brake device which is easily for moving operating table.



### Super Low Position

Super low position(500mm), suitable for operation on head.




### Kidney Bridge

Built-In Kidney Bridge, More Convenient And Nice Appearance.



### Screen



### Armrest



### Knee Curutch



### Shoulder Support






### Body Support

# Pisa 3

## Technical Chart

### Physical Parameters

Weight	1.4kg
Handle	Standard accessory, up loaded
Inter-lockable	Yes
Door	Manual
CPU	Double CPU
Ambulance use	Apply GB/T 14710-2009 transport test
IP Level	IP34







3403 NW 82nd Avenue, Suite 102 E.  
Doral, Florida 33122. USA  
[www.hosmed-inc.com](http://www.hosmed-inc.com)