



# 12-Channel Electrocardiograph **Pulsatio 12**



# 12-Channel Electrocardiograph

## Pulsatio 12


### Features

- 8.4 inch LED touch screen
- Multi-language operation and information input
- 12-lead acquisition and amplification, 12 ECG waveforms recording simultaneously
- 1:1 ratio on-screen diagnose
- Manual/ Auto/ Rhythm working modes selectable
- AC/DC power supply, built-in battery, support 4 hours working time
- 300 ECG records storage and review, built-in SD card (optional) can save 40000 cases
- Support 210mm, 216mm roll and Z-folded paper
- Net-grid on screen for easier observation

### Interface



System setup



Printing Preview




Lead status map



IP keyboard interface



Patient information interface



Review interface

### Various printing modes



12-Channel Electrocardiograph

# Pulsatio 12



8.4 inches LED touch screen



Alphanumeric Keyboard




Hidden Handle



Memory 300 ECG records;  
USB for data transmission;  
Support PC-ECG software;  
Support external printer  
connection






- Support USB printer, A4 paper printer without computer, long-term preservation and saving cost.
- Support USB flash disk back-up and recovery of patient data, integral storage expand; Support connection to computer for preserving and reading patient data.
- Support preservation and copy of cases in picture format.



12-Channel Electrocardiograph

# Pulsatio 12







3403 NW 82nd Avenue, Suite 102 E.  
Doral, Florida 33122. USA