



Infusion Pump **Frey 4.3i**

Because every drop matters,
rust Hosmed Inc.



Call Us
+1 786 266-7795



Visit Our Website
hosmed-inc.com

Safety

High accuracy: $\pm 5\%$
air bubble detect & pressure
detect Anti-bolus function
Water proof level: IP34
9 hours battery back-up time

Friendliness

4.3" colour touch screen
Intuitive operation system
Automatic calibration Inter-
stackable design



Precision in Every Drop – Infusion Pump by Hosmed Inc

Hosmed Inc. Infusion Pump integrates touch screen technology to make the operation easier and faster. With modular design, the infusion pump and injection pump can be flexibly stacked and combined into a multi-channel infusion pump, so that department resources can be adjusted reasonably at any time, and the utilization rate can be improved.

Ensure safe, accurate, and reliable delivery of fluids and medications with our state-of-the-art Infusion Pump. Designed for efficiency and patient comfort, this device combines modern technology with intuitive operation for healthcare professionals.

Call us for appointment and consultation +1 786 266-7795

Technical Specifications

Physical Parameters

Weight	1.4kg
Handle	Standard accessory, up loaded
Inter-lockable	Yes
Door	Manual
CPU	Double CPU
Pressure sensor	Up & down pressure sensors
Air bubble detector	Up & down air bubble detector (optional)
Ambulance use	Apply GB/T 14710-2009 transport test
IP Level	IP34

Display

Screen	4.3 " Color TFT touch screen
Language	English
Input interface	Touch screen Numeric & Alphanumeric keyboard & keypa

IV SET

Open/closed system	Open system/close system
IV SET	Apply to all IV sets with GB 8368 standard, with multiple brand names stored

Infusion Parameter

Infusion Mode	Rate, Time, Body-weight, drip, Drug library, Ramp up/down, Loading-dose, Sequence mode
Relay Mode	Yes
Micro mode	Yes
Infusion theory	Linear peristalsis
Rate Range	0.1-1500ml/h
Rate increment	0.01ml/h(0.01-99.99ml/h) 0.1ml/h(100-999.9ml/h) 1ml/h(1000-1500ml/h)
Infusion Accuracy	±5%
Change rate during	Yes
KVO	0.1-5ml/h adjustable
Bolus	Yes
Door	Manual/Automatic
Purge	1500ml/h

Safety

Occlusion units	mmHg/kPa/bar/psi
Occlusion levels	12 levels, 75-900 mmHg
Occlusion display	Yes
Anti-bolus	Yes
Air bubble detection	7 levels adjustable: 50ul,100ul,200ul,300ul,450ul, 600ul,800ul, default: 50ul
Anti free flow	Mechanical

Alarm

Alarm display	Audible, visible alarm with text and image
Alarm volume	10 levels adjustable
Reminder	yes

Information Database

Database capacity	More than 5000 history records, 30 last therapies
Drug library	Editable, no less than 2000 drug information
DERS	Yes

Battery

Battery Type	Rechargeable Lithium ion polymer battery
Battery life	Standard: >4.5 hours @ 25ml/h Optional: >9 hours @ 25ml/h
Charging time	≤5 hrs

Others

Automatic keypad lock	Yes
Night mode	Yes
DC Supply	Yes, 12V
Wireless LAN, IrDA	IrDA standard, Wireless LAN optional



Infusion Pump **Frey 4.3i**

Because every drop matters, trust Hosmed Inc.



Call Us
+1 786 266-7795



Visit Our Website
hosmed-inc.com