



Hydraulic Stretcher **NS5H**

Hydraulic Power. Effortless Patient Care



Call Us
+1 786 266-7795



Visit Our Website
hosmed-inc.com



Built for Safety. Designed for Performance



- Four double-sided 6-inch casters, all lockable simultaneously using the foot pedal on both sides. TPR wheels resist wear even after 30 km of use. Features an anti-winding hard shell, integrated one-piece construction with no bolts.



- Leg section lowers to raise the head section (higher than the leg section).
- Back section lowers to raise the leg section (higher than the head section).
- Step on the center of the pedal to lower the bed base smoothly and steadily.



- Waterproof PU mattress, convenient for transfer patient

Call us for appointment and consultation +1 786 266-7795

Technical Chart

QUALITY CONTROL

METAL MATERIAL

Baosteel from Fortune Global 500
Panasonic Robotics ensure 360° full
smooth welding

500 | 360°

PLASTIC

Environmental plastic can
be sterilized in 100°C,
tensile upto 30MPa

100 | 30

STEELS | PAINTING

II Process epoxy painting,
ASTM testing anti-bacterial, paint
thickness 0.12mm, brightness
60°, paint can resist 50kg impact

0.12 | 60° | 50

BRAKE

Break spring keep qualified
after 10,000times
uninteruped break

HYDRAULIC PUMPS

Height adjustment,
trendelenburg and
revised tendelenburg
through two USA
hydraulic pumps

CENTER WHEEL


Can be engaged to facilitate
direction control, enabling
users to navigate the
stretcher easily in hectic and
dynamic environments.

Technical parameters

Size	2100*850*(630-930)mm
Back-rest	0-75° (±10°)
Knee-rest	0-45°
Reverse Trendelenburg	0-15°
Trendelenburg	0-15°

Technical configuration

USA Hydraulic Pump	2 pcs.
Crank System	1 set
Back Rest Air Spring	2 pcs.
8" German Castors, Central Locking	4 pcs.
Emergency Push Handle	1 set
Folding Side Rail	1 set
74mm Wall Bumper Wheel	1 set
IV pole prevision	2 pcs.
Mattress	1 pcs.
IV Pole	1 set
Central Locking Break Pedal	1 pcs.
X-Ray Translucent Bed Platform	1 set





Hydraulic Stretcher

NS5H

Advanced Technology. Trusted Care



Call Us
+1 786 266-7795



Visit Our Website
hosmed-inc.com